Water & Bacteria Proof Film Dressing + Pad

Description

This dressing helps wounds heal with a specially designed absorbent pad which has a soft cushioning effect to protect the wound. The pad has an extra layer which minimises pain when removing the dressing. In addition, our water and bacteria-proof protective film helps reduce the risk of infection. As it is waterproof, you can shower with it. This dressing provides a mechanical barrier to water and bacteria. It does not contain bactericidal agents.

Application of dressing:

1. Wash hands before applying the dressing.
2. Clean and dry the wound and surrounding skin gently.
3. Remove dressing from pouch.
4. Remove first half of the white backing paper.
5. Apply the dressing so the white pad sits directly over the wound.
6. Remove the second white backing paper and press lightly down to secure the film to the skin.
7. Remove the top carrier where the orange arrow points.

Change the dressing when the film is not completely sticking to the skin, or you see blood/fluid seeps through the pad.

Removing the dressing:

As this is an adhesive product, carefully apply and remove the dressing especially if skin is sensitive. Gently lift one corner and stretch the clear film outwards from the pad. Continue this around the pad until dressing is removed.

Precautions for use

Seek medical advice immediately for all serious wounds and burns or if you experience redness or discomfort near the wound. Performance of product may vary if it is submerged under water for long periods.

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Single use. Store in a cool, dry place (<35 degrees Celsius).

Intended Use:

To absorb exudate and reduce the risk of infection at the applied wound area for minor wounds by providing a mechanical barrier to water and bacteria.

Basic UDI-DI - 888501891020304BASZ

Manufacturer:

BPositive Pte Ltd
www.bpositivehealthcare.com

Importer:

BPositive Pte Ltd
Singapore Business Federation Centre
160 Robinson Road, #14-04
068914 Singapore

BPositive Healthcare Pty Ltd

Suite 18, 11 Wilson Street, South Yarra VIC,
3141 Australia

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(Ver01)

Water & Bacteria Proof Film Dressing

Description

This transparent dressing helps wounds heal by acting as a protective barrier to germs and contaminants. In doing so, it helps reduce the risk of infection at the wound site. It is ideal for minor wounds like burns, cuts, grazes and covering sutures. As it is very conformable, you can use it for difficult areas like fingers, toes or even on the face. If your wound is still releasing fluid, you will need to apply an absorbent dressing directly over the wound and then cover it with this dressing. This dressing provides a mechanical barrier to water and bacteria. It does not contain bactericidal agents.

Application of dressing:

1. Wash hands before applying the dressing.
2. Clean and dry the wound and surrounding skin gently.
3. Remove dressing from pouch.
4. Remove first half of the white backing paper.
5. Apply the dressing directly over the wound.
6. Remove the top paper carrier frame leaving just the film on your skin.

Change the dressing when the film is not completely sticking to the skin.

Removing the dressing:

7. As this is an adhesive product, carefully apply and remove the dressing, especially if your skin is sensitive. Gently lift one corner and stretch the clear film outwards. Continue stretching around the dressing until it is removed.

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Intended Use:

To protect and reduce the risk of infection at the applied wound area for minor wounds by providing a mechanical barrier to water and bacteria.

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35°C / 95°F

Description

This transparent dressing helps protect catheter sites by acting as a protective barrier to germs and contaminants. In doing so, it helps reduce the risk of infection at the site.

Intended Use

The single use adhesive, polyurethane film dressing is a water & bacteria (non-antibacterial) proof film dressing, used to protect and reduce the risk of infection at the applied wound area for minor wounds, such as cuts, grazes and stitches. It may also be used to secure other devices such as vascular catheters, infusion ports or tubing.

Application of dressing:

1. Wash hands before applying the dressing.
2. Clean and dry the area and surrounding skin gently.
Use of ointments, creams, or lotions may significantly reduce adhesion to skin surface.
3. Remove dressing from pouch and apply dressing to catheter site

Removing the dressing:

4. As this is an adhesive product, carefully remove the dressing, especially if your skin is sensitive. Gently lift one corner and stretch the clear film outwards.
Continue stretching around the dressing until it is removed.

Precautions for use

Stop any bleeding at the site before applying the dressing.
Do not stretch dressing during application as tension can cause skin injury and trauma.
Allow the skin to dry thoroughly before applying the dressing to prevent skin irritation and to ensure good adhesion. Lotions and soap residue can negatively impact the adhesive and dressing performance.

Sterile, Single use: Do not reuse, reprocess or sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

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(Ver01)            35°C / 95°F

Manufacturer:
BPositive Pte Ltd (Singapore)
www.bpositivehealthcare.com

Importer

BPositive Pte Ltd (Singapore)
Singapore Business Federation Centre
160 Robinson Road, #14-04
068914 Singapore

BPositive Pty Ltd (Australia)
Suite 18, 11 Wilson Street
South Yarra VIC, 3141 Australia

T T MEDICAL MANAGEMENT SDN. BHD.
A-3A-03 EkoCheras,
No:693, Batu 5 Jalan Cheras,
56000 Kuala Lumpur, Malaysia

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Basic UDI-DI - 888501891020406BATN

Description

This transparent dressing helps protect catheter sites by acting as a protective barrier to germs and contaminants. In doing so, it helps reduce the risk of infection at the site. The dressing is bordered, notched, and reinforced with soft cloth tape to provide a better seal around catheters

Intended Use

The single use adhesive, polyurethane film dressing is a water & bacteria (non-antibacterial) proof film dressing, used to protect and reduce the risk of infection at the applied wound area for minor wounds, such as cuts, grazes and stitches. It may also be used to secure other devices such as vascular catheters, infusion ports or tubing

Application of dressing:

1. Wash hands before applying the dressing.
2. Clean and dry the area and surrounding skin gently. Use of ointments, creams, or lotions may significantly reduce adhesion to skin surface
3. Remove dressing from pouch and apply dressing to catheter site
4. After dressing has been applied, apply the sterile tape strip(s) to further secure I.V. tubing or to stabilize catheter. Refer to figures on packaging

Removing the dressing:

5. As this is an adhesive product, carefully remove the dressing, especially if your skin is sensitive. Gently lift one corner and stretch the clear film outwards. Continue stretching around the dressing until it is removed.

Precautions for use

Stop any bleeding at the site before applying the dressing.

Do not stretch dressing during application as tension can cause skin injury and trauma.

Allowed the skin to dry thoroughly before applying the dressing to prevent skin irritation and to ensure good adhesion. Lotions and soap residue can negatively impact the adhesive and dressing performance

Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.

Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Dressing is sterile unless the individual pouch is opened or damaged. Do not use the dressing if the pouch is open or damaged.

Store in a cool, dry place (<35 degrees Celsius).

Basic UDI-DI - 888501891020406BATN

Manufacturer:
BPositive Pty Ltd (Singapore)
www.bpositivehealthcare.com

Importer
BPositive Pty Ltd (Singapore)
Singapore Business Federation Centre
160 Robinson Road, #14-04
068914 Singapore

BPositive Pty Ltd (Australia)
Suite 18, 11 Wilson Street
South Yarra VIC, 3141 Australia

T T MEDICAL MANAGEMENT SDN. BHD.
A-3A-03 EkoCheras,
No:693, Batu 5 Jalan Cheras,
56000 Kuala Lumpur, Malaysia

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